



Implant Torque Specification Guide

Manufacturer	Platform	Torque (N-cm)	Driver
Astra	Astra 3.0	15	0.050" (1.27 mm) Hex
	Astra 3.5-4.0	20	0.050" (1.27 mm) Hex
	Astra 4.5-5.0	25	0.050" (1.27 mm) Hex
Biohorizons	Internal 3.5, 4.5, 5.7	30	0.050" (1.27 mm) Hex
	Tapered Internal 3.5, 4.5, 5.7	30	0.050" (1.27 mm) Hex
	External Implant 3.5	30	0.050" (1.27 mm) Hex
	External Implant 4.0, 5.0	30	0.050" (1.27 mm) Hex
Biomet 3i	Certain™ 3.25, 4.0, 4/3, 5.0, 5/4, 6.0	20	0.048" (1.22 mm) Hex
	Certain™ Prevail™ 3/4/3 (3.4), 4/5/4, 5/6/5, 4/3, 5/4, 6/5	20	0.048" (1.22 mm) Hex
	Certain™ XP 4/5, 5/6	20	0.048" (1.22 mm) Hex
	External Hex	35	0.048" (1.22mm) Hex
Blue Sky Bio	One Stage	30	Star & 0.048" Hex
	Trilobe	30	Unigrip & 0.048" Hex
	Internal Hex & Molar	30	0.050" (1.27mm) Hex
	Conus 12 & Three	30	0.050" (1.27mm) Hex
	Max	30	0.048" (1.22mm) Hex
	Quatro	30	0.048" (1.22mm) Hex
Dentsply	ANKYLOS® C/X – C/3.5, 4.5, 5.5, 7.0 (non-indexed)	15	1.0 mm Hex
	ANKYLOS® C/X – /X 3.5, 4.5, 5.5, 7.0 (indexed)	15	1.0 mm Hex
	Friadent	24	0.048" (1.22mm) Hex
	OsseoSpeed™ EV 3.0 – Green	25	0.050" Hex
	OsseoSpeed™ EV 3.6 – Purple	25	0.050" Hex
	OsseoSpeed™ EV 4.2 – Yellow	25	0.050" Hex
	OsseoSpeed™ EV 4.8 – Blue	25	0.050" Hex
	OsseoSpeed™ EV 5.4 – Brown	25	0.050" Hex
	OsseoSpeed™ Profile EV 4.2 Yellow	25	0.050" Hex
	OsseoSpeed™ Profile EV 4.8 Blue	25	0.050" Hex
	OsseoSpeed™/OsseoSpeed™ TX 3.0 – Yellow	15	0.050" Hex
	OsseoSpeed™/OsseoSpeed™ TX 3.5, 4.0 – Aqua	20	0.050" Hex
	OsseoSpeed™/OsseoSpeed™ TX 4.5, 5.0 – Lilac	25	0.050" Hex
	OsseoSpeed™ 4.0 – Blue (Previous)	25	0.050" Hex
	OsseoSpeed™ TX Profile 4.5, 5.0	25	0.050" Hex
	XIVE® S 3.0	24	1.22 mm Hex
	XIVE® S 3.4, 3.8, 4.5, 5.5	24	1.22 mm Hex
FRIALIT™ 3.4, 3.8, 4.5, 5.5	24	1.22 mm Hex	
Glidewell/Hann	iNCLUSIVE® Tapered 3.7, 4.7, 5.2	35	0.050" (1.27mm) Hex
Neodent	CM 15	15	0.035" (0.9mm) Hex
	CM 32	32	0.048" (1.22mm) Hex
	CM 32	32	0.063" (1.6mm) Hex
Nobel Biocare	Angulated Screw Channel	35	Nobel Biocare® Omnigrip
	NobelActive™ 3.0	15	Nobel Biocare® Unigrip
	NobelActive™/NobelReplace™ Conical Connection NP 3.5	35	Nobel Biocare® Unigrip
	NobelActive™/NobelReplace™ Conical Connection RP 4.3, 5.0	35	Nobel Biocare® Unigrip
	NobelReplace™ NP 3.5	35	Nobel Biocare® Unigrip
	NobelReplace™ RP 4.3	35	Nobel Biocare® Unigrip
	NobelReplace™ WP 5.0	35	Nobel Biocare® Unigrip
	NobelReplace™ 6.0	35	Nobel Biocare® Unigrip
	Brånemark System™ NP 3.3	35	Nobel Biocare® Unigrip
	Brånemark System™ RP 3.75, 4.0	35	Nobel Biocare® Unigrip
	Brånemark System™ WP 5.0, 5.5	35	Nobel Biocare® Unigrip
Straumann	Bone Level 3.3 NC (Narrow CrossFit™)	35	SCS Screwdriver
	Bone Level 4.1, 4.8 RC (Regular CrossFit™)	35	SCS Screwdriver
	Standard/Standard Plus 4.8 RN (Regular Neck)	35	SCS Screwdriver
	Standard/Standard Plus 4.8 WN (Wide Neck)	35	SCS Screwdriver
	Standard Plus 3.3 NN (Narrow Neck)	35	SCS Screwdriver
	Tissue Level	30	Star/Torx
Zimmer/Zimvie	Tapered Screw-Vent 3.5, 4.5, 5.7/Tapered Screw-Vent 3.5, 4.5, 5.7 friction fit	30	(.050) 1.25 mm Hex
	Screw-Vent 3.3, 3.7, 4.5/Screw-Vent 3.3, 3.7, 4.5 friction fit	35	1.25 mm Hex
	SwissPlus™ 4.8	35	1.25 mm Hex